

# PFAS Analysis by EPA 1633 | Building a Comprehensive Capability

**Waters Corporation**

July 10, 2024



## Today's presenters



**Stuart Oehrle**

Senior Principal Technical  
Support Specialist



**Chelsea Plummer**

Senior Product Marketing  
Manager



**Jennifer Burzynski**

Field Service Applications  
Team Leader

# Waters at a glance

- Liquid chromatography, mass spectrometry, and thermal analysis innovator
- Industry-leading service and informatics
- Partner with customers around the world to ensure the:
  - Purity of our water and food
  - Efficacy, safety & quality of medicines and vaccines
  - Durability of products we use everyday



# \$2.9BN

2023 Revenue

# ~7,900

Employees Worldwide

# 35+

Countries with Operations



MEMBER OF  
**Dow Jones Sustainability Indices**  
Powered by the S&P Global CSA



# What we will cover today

- Getting started with PFAS testing: what you need for comprehensive testing capability
  - Spotlight on LC-MS/MS technology
  - Spotlight on sample preparation
- Training and support as foundations for success
- Q&A




# What we will cover today

- Getting started with PFAS testing: what you need for comprehensive testing capability
  - Spotlight on LC-MS/MS technology
  - Spotlight on sample preparation
- Training and support as foundations for success
- Q&A



# Reminder of why (and what) we are testing: EPA method 1633






Newly finalized method (January 2024) to test for **40 PFAS compounds** in **wastewater, surface water, groundwater, soil, biosolids, sediment, landfill leachate, and fish tissue**



Created by the **US EPA's Office of Water**, partnership with the **US Department of Defense's Strategic Environmental Research and Development Program (SERDP)**

Validated through a **multi-laboratory study** to finalize method and add performance criteria


Intended to be used for:

-  Clean Water Act Compliance (wastewater discharge permits – NPDES)
-  Superfund sites
-  Remediation and investigation programs (including for DoD installations)

# Brief intro to liquid chromatography + tandem mass spectrometry (LC-MS/MS): the LC




# Brief intro to liquid chromatography + tandem mass spectrometry (LC-MS/MS): the MS







# LC-MS/MS has been selected by regulators for PFAS for the technology's selectivity and sensitivity

- High selectivity ➤ Reduce or eliminate matrix interferences
- High sensitivity ➤ Trace levels - low reporting limits
- Quantitative accuracy ➤ Reproducibility, stability, dynamic range and accurate quantitation of targets at low levels in matrix
- Robustness ➤ Complex sample matrices, reduced sample clean-up




Waters ACQUITY™ Premier  
UPLC™ + Xevo™ TQ Absolute MS

# What does that give you?



# The full picture




**46 native PFAS** covering both EPA 1633 and ASTM 8421 lists

# PFAS analysis by EPA 1633 – results from LC-MS/MS


Certified reference material (CRM) validation

## PFAS in wastewater CRM



Mean trueness of 92%  
Trueness range of 73 – 112%


## PFAS in soil CRM




Mean trueness of 97%  
Trueness range of 85 – 120%

# LC-MS/MS is also used for EPA PFAS drinking water methods

## EPA 537.1



## EPA 533



Identification and amounts (typically in ng/L)


## Putting ng/L into perspective

Nanograms of compound per liter of solution (aka part-per-trillion (ppt))

– **One** part per 1,000,000,000,000 ( $10^{12}$ ) parts

### Some ways to visualize one ng/L ....


1 drop of water in 20 Olympic-sized swimming pools



One second of time in approximately 31,700 years



Traveling 6 inches out of a 93-million-mile journey to the sun<sup>1</sup>




1. The earth is 93 million miles (avg.) from the sun

# Testing for PFAS requires a comprehensive analysis solution



## Sample prep

- SPE cartridges
- Sample preparation apparatus (manual vacuum manifold or an automated system)
- Nitrogen sample evaporator




## Instruments

- Ultra-Performance Liquid Chromatography (UPLC)
- Tandem quadrupole mass spectrometer (MS/MS)
- PFAS hardware kit
- Analytical columns



## Informatics

- Instrument control and acquisition software
- Data processing and reporting
- Output to LIMS system or similar database
- Local quality system



## CRM and PT

- Certified reference materials (CRM)
- Proficiency testing (PT) program for analysts

# What we will cover today




- Getting started with PFAS testing: what you need for comprehensive testing capability
  - Spotlight on LC-MS/MS technology
  - Spotlight on sample preparation
- Training and support as foundations for success
- Q&A





# Solid-phase extraction (SPE) - why we use it for PFAS analysis

- **Regulatory compliance**
  - US EPA 533, 537.1, and 1633 require solid phase extraction clean up steps
  - Reproducibility across laboratories running same method
- **Clean up from complex matrices**
  - Environmental waters, solids, biosolids, and other samples contain interferences
  - Reduce ion suppression and background
- **Concentration of PFAS**
  - Measuring at low levels (ng/L)
  - Increase limits of detection
- **Selectivity in extraction**
  - Use of WAX sorbent more specific to PFAS



# Sample preparation for EPA Method 1633

EPA Method 1633 requires cleanup with both a weak-anion exchange (WAX) SPE cartridge AND carbon

Depending on the sample matrix, the order of clean up is determined

EPA 1633 is a **performance-based** method, allowing for modifications as long as recoveries and RSDs% meet the method criteria



WAX SPE cartridges and loose graphitized carbon black (GCB)

**Ground, surface, and wastewaters**



WAX cleanup first, then GCB

**Soils, sediments, biosolids, tissue**



GCB cleanup first, then WAX



WAX/GCB and GCB/WAX dual-phase, or bilayer cartridges combine these two cleanup steps into one single cartridge for faster and easier use in EPA 1633



# SPE protocol for EPA 1633 using dual-phase cartridges

## 1 Prepare sample

- Some Pre-SPE steps dependent on sample matrix type
- Check pH and adjust to  $6.5 \pm 0.5\%$  using formic acid or ammonium hydroxide as needed
- Check cartridge orientation, pack glass wool to half height of the SPE barrel

## 2 Condition and equilibrate SPE cartridge

- 15 mL of 1% (v/v) ammonium hydroxide methanol
- 5 mL of 0.3 M formic acid

## 3 Load sample onto SPE cartridge and washes

- Load sample at 5 mL/min
- Wash cartridge with 10 mL of reagent water, ensuring to rinse reservoir with this solution
- Wash with 5 mL of 1:1 0.1M formic acid:methanol, ensuring to rinse the reservoir with this solution
- Dry cartridge for 15 seconds

## 4 Elute sample from SPE cartridge

- Place collection tubes in manifold
- Rinse bottle with 5 mL 1% (v/v) ammonium hydroxide in methanol. Transfer to cartridge and elute
- Add 25  $\mu\text{L}$  acetic acid to each sample
- Spike each sample with non-extracted internal standard





Full workflow with  
WAX/GCB and  
GCB/WAX  
cartridges video

# Sample prep that meets EPA Method 1633 requirements

- Proof of equivalence *must* be demonstrated when using a dual-phase cartridge.
- Dual-phase cartridges remove 4 steps in the protocol, and are faster, easier, and cleaner to use than loose GCB.
- **Recoveries and repeatabilities** for GCB/WAX for soil and fish tissue and WAX/GCB for ground, surface, influent and effluent waters **were measured and are within acceptance criteria and are repeatable.**



- For soil (n=3):
  - Mean recovery of all EIS was 81%
  - Mean %RSD was 2.8%
- For tissue (n=3):
  - Mean recovery of all EIS was 85%
  - Mean %RSD was 9.2%
- For environmental waters (n=5):
  - Mean recovery of all EIS was 91.2%
  - Mean %RSD was 9.2%



# Looking ahead - Making analysis by EPA 1633 'easier'

## Research underway

---

Fully automate SPE using Oasis WAX/GCB using PromoChrom SPE-03

Introduction of dual-phase cartridges with WAX and GCB to reduce manual carbon steps

Utilize the high sensitivity of the Xevo TQ Absolute to reduce aqueous sample sizes



PromoChrom  
SPE-03



Oasis WAX/GCB  
and GCB/WAX







ACQUITY Premier with  
Xevo TQ Absolute

# What we will cover today

- Getting started with PFAS testing: what you need for comprehensive testing capability
  - Spotlight on LC-MS/MS technology
  - Spotlight on sample preparation
- Training and support as foundations for success
- Q&A



# Decoding PFAS literature



NMeFOSAA

LFSM

SUR

<sup>13</sup>C-PFHxA

EIS

NIS

CCC

ADONA





## What do the methods *not* give us?

- Instrument-specific information
  - Methods
  - Parameters
- Software-specific information
  - Data processing
  - Calculations
- A partner



# Education and training are fundamental to successfully establishing a PFAS analysis capability

Waters™

## Example suite of training options

- Custom onsite training
- Virtual training
- Advanced method consultation
- Education subscriptions
- Certified in-house courses
- E-learning courses




**No matter what you choose, we support you.**

The Waters™ Customer Education Team is here for you with expert instructors whether you have just purchased a new Waters instrument, need help learning a new software or application, or simply require refresher training on any of our platforms.

Support for PFAS analysis should span from preparing your lab through the ability to run the analysis independently

Waters™




We are with you every step of the way

## PFAS Solution Setup Guide

- Consultant-led guidance
- Ensures necessary preparation for PFAS analysis

## PFAS Application Success Guide

- Detailed reference document
- Instructions for entire end-to-end workflow



60 mL Sample Reservoir  
Sample Reservoir Adaptor  
SPE Cartridge  
SPE Manifold  
Waste

Z1887723 V01R 00 32 OF 41

Sample	Isomer	Concentration (ng/mL)	Report Value (ng/mL)
Standard 1	Perfluorooctanoic acid	1000	1000
Standard 2	Perfluorooctanoic acid	2000	2000
Standard 3	Perfluorooctanoic acid	5000	5000
Standard 4	Perfluorooctanoic acid	10000	10000
Standard 5	Perfluorooctanoic acid	20000	20000
Standard 6	Perfluorooctanoic acid	50000	50000
Standard 7	Perfluorooctanoic acid	100000	100000
Standard 8	Perfluorooctanoic acid	200000	200000
Standard 9	Perfluorooctanoic acid	500000	500000
Standard 10	Perfluorooctanoic acid	1000000	1000000
Standard 11	Perfluorooctanoic acid	2000000	2000000
Standard 12	Perfluorooctanoic acid	5000000	5000000
Standard 13	Perfluorooctanoic acid	10000000	10000000
Standard 14	Perfluorooctanoic acid	20000000	20000000
Standard 15	Perfluorooctanoic acid	50000000	50000000
Standard 16	Perfluorooctanoic acid	100000000	100000000
Standard 17	Perfluorooctanoic acid	200000000	200000000
Standard 18	Perfluorooctanoic acid	500000000	500000000
Standard 19	Perfluorooctanoic acid	1000000000	1000000000
Standard 20	Perfluorooctanoic acid	2000000000	2000000000
Standard 21	Perfluorooctanoic acid	5000000000	5000000000
Standard 22	Perfluorooctanoic acid	10000000000	10000000000
Standard 23	Perfluorooctanoic acid	20000000000	20000000000
Standard 24	Perfluorooctanoic acid	50000000000	50000000000
Standard 25	Perfluorooctanoic acid	100000000000	100000000000
Standard 26	Perfluorooctanoic acid	200000000000	200000000000
Standard 27	Perfluorooctanoic acid	500000000000	500000000000
Standard 28	Perfluorooctanoic acid	1000000000000	1000000000000
Standard 29	Perfluorooctanoic acid	2000000000000	2000000000000
Standard 30	Perfluorooctanoic acid	5000000000000	5000000000000
Standard 31	Perfluorooctanoic acid	10000000000000	10000000000000
Standard 32	Perfluorooctanoic acid	20000000000000	20000000000000
Standard 33	Perfluorooctanoic acid	50000000000000	50000000000000
Standard 34	Perfluorooctanoic acid	100000000000000	100000000000000
Standard 35	Perfluorooctanoic acid	200000000000000	200000000000000
Standard 36	Perfluorooctanoic acid	500000000000000	500000000000000
Standard 37	Perfluorooctanoic acid	1000000000000000	1000000000000000
Standard 38	Perfluorooctanoic acid	2000000000000000	2000000000000000
Standard 39	Perfluorooctanoic acid	5000000000000000	5000000000000000
Standard 40	Perfluorooctanoic acid	10000000000000000	10000000000000000
Standard 41	Perfluorooctanoic acid	20000000000000000	20000000000000000
Standard 42	Perfluorooctanoic acid	50000000000000000	50000000000000000
Standard 43	Perfluorooctanoic acid	100000000000000000	100000000000000000
Standard 44	Perfluorooctanoic acid	200000000000000000	200000000000000000
Standard 45	Perfluorooctanoic acid	500000000000000000	500000000000000000
Standard 46	Perfluorooctanoic acid	1000000000000000000	1000000000000000000
Standard 47	Perfluorooctanoic acid	2000000000000000000	2000000000000000000
Standard 48	Perfluorooctanoic acid	5000000000000000000	5000000000000000000
Standard 49	Perfluorooctanoic acid	10000000000000000000	10000000000000000000
Standard 50	Perfluorooctanoic acid	20000000000000000000	20000000000000000000
Standard 51	Perfluorooctanoic acid	50000000000000000000	50000000000000000000
Standard 52	Perfluorooctanoic acid	100000000000000000000	100000000000000000000
Standard 53	Perfluorooctanoic acid	200000000000000000000	200000000000000000000
Standard 54	Perfluorooctanoic acid	500000000000000000000	500000000000000000000
Standard 55	Perfluorooctanoic acid	1000000000000000000000	1000000000000000000000
Standard 56	Perfluorooctanoic acid	2000000000000000000000	2000000000000000000000
Standard 57	Perfluorooctanoic acid	5000000000000000000000	5000000000000000000000
Standard 58	Perfluorooctanoic acid	10000000000000000000000	10000000000000000000000
Standard 59	Perfluorooctanoic acid	20000000000000000000000	20000000000000000000000
Standard 60	Perfluorooctanoic acid	50000000000000000000000	50000000000000000000000
Standard 61	Perfluorooctanoic acid	100000000000000000000000	100000000000000000000000
Standard 62	Perfluorooctanoic acid	200000000000000000000000	200000000000000000000000
Standard 63	Perfluorooctanoic acid	500000000000000000000000	500000000000000000000000
Standard 64	Perfluorooctanoic acid	1000000000000000000000000	1000000000000000000000000
Standard 65	Perfluorooctanoic acid	2000000000000000000000000	2000000000000000000000000
Standard 66	Perfluorooctanoic acid	5000000000000000000000000	5000000000000000000000000
Standard 67	Perfluorooctanoic acid	10000000000000000000000000	10000000000000000000000000
Standard 68	Perfluorooctanoic acid	20000000000000000000000000	20000000000000000000000000
Standard 69	Perfluorooctanoic acid	50000000000000000000000000	50000000000000000000000000
Standard 70	Perfluorooctanoic acid	100000000000000000000000000	100000000000000000000000000
Standard 71	Perfluorooctanoic acid	200000000000000000000000000	200000000000000000000000000
Standard 72	Perfluorooctanoic acid	500000000000000000000000000	500000000000000000000000000
Standard 73	Perfluorooctanoic acid	1000000000000000000000000000	1000000000000000000000000000
Standard 74	Perfluorooctanoic acid	2000000000000000000000000000	2000000000000000000000000000
Standard 75	Perfluorooctanoic acid	5000000000000000000000000000	5000000000000000000000000000
Standard 76	Perfluorooctanoic acid	10000000000000000000000000000	10000000000000000000000000000
Standard 77	Perfluorooctanoic acid	20000000000000000000000000000	20000000000000000000000000000
Standard 78	Perfluorooctanoic acid	50000000000000000000000000000	50000000000000000000000000000
Standard 79	Perfluorooctanoic acid	100000000000000000000000000000	100000000000000000000000000000
Standard 80	Perfluorooctanoic acid	200000000000000000000000000000	200000000000000000000000000000
Standard 81	Perfluorooctanoic acid	500000000000000000000000000000	500000000000000000000000000000
Standard 82	Perfluorooctanoic acid	1000000000000000000000000000000	1000000000000000000000000000000
Standard 83	Perfluorooctanoic acid	2000000000000000000000000000000	2000000000000000000000000000000
Standard 84	Perfluorooctanoic acid	5000000000000000000000000000000	5000000000000000000000000000000
Standard 85	Perfluorooctanoic acid	10000000000000000000000000000000	10000000000000000000000000000000
Standard 86	Perfluorooctanoic acid	20000000000000000000000000000000	20000000000000000000000000000000
Standard 87	Perfluorooctanoic acid	50000000000000000000000000000000	50000000000000000000000000000000
Standard 88	Perfluorooctanoic acid	100000000000000000000000000000000	100000000000000000000000000000000
Standard 89	Perfluorooctanoic acid	200000000000000000000000000000000	200000000000000000000000000000000
Standard 90	Perfluorooctanoic acid	500000000000000000000000000000000	500000000000000000000000000000000
Standard 91	Perfluorooctanoic acid	1000000000000000000000000000000000	1000000000000000000000000000000000
Standard 92	Perfluorooctanoic acid	2000000000000000000000000000000000	2000000000000000000000000000000000
Standard 93	Perfluorooctanoic acid	5000000000000000000000000000000000	5000000000000000000000000000000000
Standard 94	Perfluorooctanoic acid	10000000000000000000000000000000000	10000000000000000000000000000000000
Standard 95	Perfluorooctanoic acid	20000000000000000000000000000000000	20000000000000000000000000000000000
Standard 96	Perfluorooctanoic acid	50000000000000000000000000000000000	50000000000000000000000000000000000
Standard 97	Perfluorooctanoic acid	100000000000000000000000000000000000	100000000000000000000000000000000000
Standard 98	Perfluorooctanoic acid	200000000000000000000000000000000000	200000000000000000000000000000000000
Standard 99	Perfluorooctanoic acid	500000000000000000000000000000000000	500000000000000000000000000000000000
Standard 100	Perfluorooctanoic acid	1000000000000000000000000000000000000	1000000000000000000000000000000000000

These are pre-calculated QC samples. Adding known concentration value here in ng/mL will automatically calculate percent recovery. User Divisor for these are set to 1.0

Insert calculated User Factor here before processing to adjust for determined sample volume collected.

User Divisor of 0.25 report values in ng/L to 250 to report as ng/mL if desired.

Figure 8 - Example sample list with parameters assigned for data processing. User Factor is calculated in the accompanying User Guide Excel file (8) to account for the exact sample volume collected. Using a User Divisor of 250 will report values in ng/mL units.


### 13.2 Stock concentration factors

Some PFAS compounds (usually those with sulfonate groups) are present in standard solutions in their salt form. These compounds are detected on the mass spectrometer in their acid form. Therefore, the standard manufacturer will typically provide two certified concentration levels on the certificate of analysis (COA) to provide the customers with concentrations of these compounds in both the salt and the acid forms.

Also, some PFAS compounds are found in standards in both their linear and branched isomer forms. This changes the overall concentration of each compound in the standard mix that has isomers present. The nominal concentration listed on the standard vial (for example, 2000 ng/mL) is the total concentration of all isomers. To properly calculate the concentration of each isomer group separately, the ratio of each isomer (provided in the COA, example shown in Figure 9) must be considered.

Both the isomer concentration and the acid concentration can be automatically taken into account during the sample processing when using TargetLynx. By using the Stock Concentration Factor parameter in the TargetLynx method editor, the Stock Concentration Factor is a percentage and will automatically apply the designated percentage against the specified concentration in the sample list (Conc A, Conc B, and so on).

# PFAS application user training: recap



# What we will cover today

- Getting started with PFAS testing: what you need for comprehensive testing capability
  - Spotlight on LC-MS/MS technology
  - Spotlight on sample preparation
- Training and support as foundations for success
- Q&A




What questions  
do you have?



# Want to learn more? Additional resources


## Regulatory method and application notes



**EPA** United States Environmental Protection Agency  
Office of Water  
www.epa.gov January 2024

**Method 1633**  
**Analysis of Per- and Polyfluoroalkyl Substances (PFAS) in Aqueous, Solid, Biosolids, and Tissue Samples by LC-MS/MS**


[EPA Method 1633](#)



**Application Note**  
Analysis of Per- and Polyfluoroalkyl Substances (PFAS) in Accordance With EPA 1633 Part 1: Establishing and Assessing the Method

Karl L. Organtini, Kenneth J. Rosnack, Peter Hancock  
Waters Corporation

[Waters EPA 1633 Application Note Part 1](#)



**Application Note**  
Analysis of Per- and Polyfluoroalkyl Substances (PFAS) in Accordance With EPA 1633 Part 2: Analysis of Aqueous Matrices

Karl L. Organtini, Kenneth J. Rosnack, Chelsea Plummer, Peter Hancock, Oliver Burt  
Waters Corporation

[Waters EPA 1633 Application Note Part 2](#)



[ WHITE PAPER ]


Waters™

**Best practices for monitoring PFAS contamination in a routine shared-space commercial laboratory**

Nicola Dreelin, Henry Fobdy, Karl Organtini, Stuart Adams, Karl Rosnack, Peter Hancock  
Waters Corporation, Milford, MA, USA and Waters Corporation, Winslow, UK

[Best practices for monitoring PFAS contamination](#)

## Clean Water State Revolving Fund (CWSRF) resources



**What Types of Clean Water Projects are Eligible?**

Planning, design, and construction activities that are eligible for the base CWSRF and CWASA program, specifically address one or more ECS, and result in a water quality benefit are eligible. Project types include:

- ✓ Laboratory and sampling equipment purchase
- ✓ Wastewater treatment (centralized & decentralized)
- ✓ Biosolids treatment/disposal
- ✓ Septic-to-sewer conversions
- ✓ Water reuse
- ✓ Source water protection
- ✓ Lagoon improvements and decommissioning
- ✓ Nonpoint source pollution control
- ✓ Stormwater management and treatment
- ✓ Landfill capping and leachate control
- ✓ Cleanup of contaminated sites
- ✓ Groundwater and surface water protection and restoration

[CWSRF Emerging Contaminants Home Page](#)

**State CWSRF Program Contacts**

CWSRF assistance is provided directly from state agencies. Contact the CWSRF program in your state for information on how to apply.

State	Contact	Phone	Email	State CWSRF Website
Alabama	Eric Rensly	(334) 271-7805	<a href="mailto:eric_rensly@adm.alabama.gov">eric_rensly@adm.alabama.gov</a>	<a href="#">Alabama Department of Environmental Management</a>
Alaska	Carrie Bohan	(907) 465-5143	<a href="mailto:carrie.bohan@alaska.gov">carrie.bohan@alaska.gov</a>	<a href="#">Alaska Department of Environmental Conservation</a>
	Young Ha	(907) 269-7544	<a href="mailto:young.ha@alaska.gov">young.ha@alaska.gov</a>	
Arizona	Lindsay Jones	(602) 364-1324	<a href="mailto:ljjones@azwfa.gov">ljjones@azwfa.gov</a>	<a href="#">Water Infrastructure Finance Authority of Arizona</a>
	General	(602) 364-1310	<a href="mailto:contact@azwfa.gov">contact@azwfa.gov</a>	
Arkansas	Debbie Dickson	(501) 682-6546	<a href="mailto:debbie.dickson@arkansas.gov">debbie.dickson@arkansas.gov</a>	<a href="#">Arkansas Department of Agriculture</a>

[State CWSRF websites and contacts](#)




Source: US Environmental Protection Agency, Waters Corporation



# Conquer your analytical PFAS challenges today

- Visit our [PFAS Resource Hub](#) for a complete guide to PFAS analysis

## PFAS Application Notebook



[DOWNLOAD NOW](#)

## Case Study: Lower Regulatory Limits Prompt Innovation in PFAS Testing with SGS



[DOWNLOAD NOW](#)

## Infographic: Oasis SPE for PFAS Testing Saves Time and Reduces Risk



[DOWNLOAD NOW](#)



Sample Collection



Sample Preparation



Chromatographic Separation



Mass Spectrometry



Data Analysis



Professional Services



ERA PT & CRM